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Arterial pressure

Arterial pressure is defined as the hydrostatic pressure exerted by the blood over the arteries as a result of the heart left ventricle contraction. Systolic arterial pressure is the higher blood pressure reached by the arteries during systole (ventricular contraction), and diastolic arterial pressure is the lowest blood pressure reached during diastole (ventricular relaxation). In a healthy young adult at rest, systolic arterial pressure is around 110 mmHg and diastolic arterial pressure is around 70 mmHg.

Blood flow is the blood volume that flows through any tissue in a determined period of time (typically represented as ml/min) in order to bring tissue oxygen and nutrients transported in blood. Blood flow is directly affected by the blood pressure as blood flows from the area with more pressure to the area with less pressure. Greater the pressure difference, higher is the blood flow. Blood is pumped from the left ventricle of the heart out to the aorta where it reaches its higher pressure levels.

Blood pressure falls as blood moves away from the left ventricle until it reaches 0 mm Hg, when it returns to the heart's right atrium. Figure 1 represents pressure changes.




Figure 1. Pressure changes on blood vessels

Blood pressure monitor operating principle

Blood pressure monitor operation is based on the oscillometric method. This method takes advantage of the pressure pulsations taken during measurements. An occluding cuff is placed on the left arm and is connected to an air pump and a pressure sensor. Cuff is inflated until a pressure greater than the typical systolic value is reached, then the cuff is slowly deflated. As the cuff deflates, when systolic pressure value approaches, pulsations start to appear. These pulsations represent the pressure changes due to heart ventricle contraction and can be used to

calculate the heartbeat rate. Pulsations grow in amplitude until mean arterial pressure (MAP) is reached, then decrease until they disappear.

Oscillometric method determines the MAP by taking the cuff pressure when the pulse with the largest amplitude appears. Systolic and diastolic values are calculated using algorithms that vary among different medical equipment developers.

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