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Pulse Oximeters in clinics have a finger clip type probe that has an LED on one side and a photo detector on the other side. The light emitted from one side of the finger travels through tissue, venous blood and arterial blood and is collected in the detector. Most of the light is absorbed or scattered before it reaches the photo detector in the other side of the finger. The flow of blood is heartbeat induced, or pulsatile in nature so the transmitted light changes with time. Red and infrared lights are used for pulse oximetry to estimate the true hemoglobin oxygen saturation of arterial blood. Oxyhemoglobin (HbO2) absorbs visible and infrared (IR) light differently than deoxyhemoglobin (Hb), and appears bright red as opposed to the darker brown Hb. Absorption in the arterial blood is represented by an AC signal which is superimposed on a DC signal representing absorptions in other substances like pigmentation in tissue, venous, capillary, bone, and so forth. Cardiac-synchronized AC signal is approximately 1% of the DC level. This is referred to as the perfusion index %. The ratio of ratios 'R' is approximated in Equation 1. % SpO2 is calculated as follows:

$$R = (ACrms\ of\ Red / DC\ of\ Red) / (ACrms\ of\ IR / DC\ of\ IR) \quad (1)$$

The standard model of computing SpO2 is defined as shown in Equation 2. This model is often used in the literature in the context of medical devices. However, accurate % SpO2 is computed based on the empirical calibration of the ratio of ratios for the specific device.

$$\% SpO2 = 110 - 25 \times R$$


Figure 1. Red and Infra Red Modulated by Cycling Blood



Figure 2. Red/Infrared Modulation Ratio

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